

MB27



- The MB27 provides accurate results at 0.01%, 0.001g. With a large capacity of 90 g.
- Simple to operate, Just press and hold the temperature or time buttons to set drying parameters.
- The compact footprint (17 × 14 × 28 cm) of the OHAUS MB27 is designed to take up less space, and easy to clean.



Moisture Analysis



Backlight



RS232 Interface



MB25

- The MB25 features a bright backlit LCD and easy one-button setup for simple use and menu navigation.
- The MB25 utilizes a halogen heating element for swift drying of samples.
- The MB25's design makes it easy-to-clean and its compact size allows it to take up less space and fit into almost any work setting.



Moisture Analysis



Backlight



RS232 Interface

Model	MB27	MB25	MB23
Capacity (g)	90	110	110
Readability	0.01%/0.001g	0.05% / 0.005g	0.1% / 0.01g
Repeatability (g)	0.2% (3 g sample) 0.05% (10 g sample)	0.2% (3 g sample) 0.05% (10 g sample)	0.3% (3g sample) 0.2% (10g sample)
Pan Size Diameter(mm)	Ø 90		Ø 90
Timer	1 - 99 minutes, 30-second increments to 60 minutes		
Heating Element	Halogen		IR
Temperature Range	50° to 160° C in 1° C increments		
Dimensions (WxHxD) (mm)	170 x 130 x 280		
Displays	% moisture or % solids or weight (g), temperature, time		

MB Accessories	
OHAUS #	Description
11113857	Temperature Calibration Kit for all MB Models
30045641	Printer Impact SF40A
12120799	Paper reels for Printer SF40A
12120798	Ink Ribbon for Printer SF40A
80252475	In-Use Cover for MB25 and MB23
80850086	Pan, Aluminum, Set(80), MB Series
80850087	Filter, Glass Fiber, Set(200)
80850043	Anti-Theft Device, MB AX EX
80252476	Holder, Pan, MB Series, MB23/MB25/MB27/MB35/MB35
80252477	Cage, Sample, MB Series
80252478	Pan Reusable 7mm 3Pcs MB Series
80252479	Pan Reusable 14mm 3Pcs MB Series

DISTRIBUTOR :-

Blue Horizon - Gopal Niwas, 135, Princess Street, Mumbai - 400 002. INDIA • Mob.: +91-98202 06611 • Email: sales.bluehorizon@gmail.com